## CombHeliX® FSHD probes

REF FSHD-HYB-001

### Instructions for use

Version: 04

Date: November 2015

#### **Background information**

FSHD (Facio-Scapulo-Humeral Muscular Dystrophy) is a genetic disease with autosomal dominant transmission defined by a progressive weakening and loss of skeletal muscles of the face (facio), the shoulder girdle (scapulo) and the upper arms (humeral). In 95 % of FSHD patients (also referred as FSHD1), a contracted 3.3kb-D4Z4 repeat array (1 to 10 copies) is present on a chromosome 4q35.2 subtelomere [1, 2] whereas the array contains between 11 and 150 repeats in the control population [3, 4]. In contrast, the FSHD2 form which is a rare form of FSHD is not linked to D4Z4 repeat array contraction. D4Z4 repeat arrays are also found on chromosome 10q26.1 subtelomeres, with no apparent relation to the disease. Chromosomes 4q and 10q subtelomeres may bear the haplotype A or B, which show several kilobases of divergent sequences. D4Z4 contraction is only pathogenic when it segregates in *cis* with the haplotype A on chromosome 4 [5-7]. The FSHD probes enable the determination of the size of the D4Z4 repeat array and the A or B haplotype on both chromosomes 4q35.2 and 10q26.1 [8].


#### **Intended use**

The “CombHeliX® FSHD probes” are designed for the diagnostic of FSHD1 by enabling the characterization of the 3.3kb-D4Z4 repeat-containing loci on chromosomes 4 and 10 by fluorescent hybridization on combed DNA extracted from blood samples and prepared according to the Molecular Combing procedures [8]. The “CombHeliX® FSHD probes” are intended to be used for *in vitro* diagnostic purpose.

Important note: In addition to the “CombHeliX® FSHD probes”, the diagnostic of FSHD requires the mandatory use of the CE-IVD “CombHeliX® FSHD Software” for the

detection of the hybridized “CombHeliX® FSHD probes”, the analysis and the interpretation of the results (see the Instruction for use for “CombHeliX® FSHD Software”).

### Principle of the diagnostic of FSHD

Molecular Combing allows the direct visualization and analysis of single DNA molecules. It enables DNA (extracted from fresh blood cells of a patient for the *in vitro* diagnostic of FSHD: step 1) to be stretched uniformly and irreversibly attached to a treated glass coverslip (step 2) according to the protocols for “DNA extraction and Molecular Combing” that are provided by Genomic Vision. The presence of specific genomic sequences (subtelomere regions of 4q35.2 and 10q26.1 in the test material) can be detected on stretched DNA by hybridization with labeled polynucleotides (“CombHeliX® FSHD probes”: step 3) and visualized by fluorescence microscopy (step 4) as it is described in the present Instructions for Use. The constant stretching factor allows the direct measurement of the probe length and the physical cartography of the region of interest. For the *in vitro* diagnostic of FSHD, the number of D4Z4

repeats in the test material is determined with the “CombHeliX® FSHD Software” (step 5) as described in the corresponding Instructions for Use provided by Genomic Vision.

### **Product content and description**

“CombHeliX® FSHD probes” set: 5 x 20 µl per vial (2 tests per vial, 10 tests in total). One test is defined as sufficient for one hybridization experiment of combed DNA on coverslips (22 x 22 mm area).

The “CombHeliX® FSHD probes” contain fluorescein-, digoxigenin- or/and biotin-labeled polynucleotides that are intended to be used in combination with specific reagents (see Reagents and materials recommended but not provided) that allow their detection in green, blue, red and magenta fluorescent signals.

Labeled polynucleotides are premixed with blocking DNA in formamide-free hybridization buffer (Saline-Sodium Citrate (SSC), Sodium Dodecyl Sulfate (SDS), Sodium Lauroyl Sarcosinate and BlockingAid™ blocking solution).

### **Storage and handling**

Upon arrival, the product must be stored between -25°C and -10°C protected from the light until the expiry date printed on the label. Improper storage of the product can destroy or impair the performance of the product and consequently the *in vitro* diagnostic assay should not be performed from such reagent since it may affect the result of the test.

Handle all reagents and slides containing fluorophores in reduced light to prevent photobleaching.

Once thawed and prior to opening, the vials should be briefly centrifuged to ensure the contents are collected at the bottom of the vials.


Once opened, the remaining content of the vial can be frozen again and stored between -25°C and -10°C up to 6 months.

### **Warnings and precautions**

**For *in vitro* diagnostic use. For professional use only.** Carefully read the operating instructions before use.

The diagnostic performances of the “CombHeliX® FSHD probes” were validated using coverslips with very long genomic DNA molecules linearly stretched at high density (picture 5) extracted from fresh blood samples. Using coverslips with genomic DNA at

lower density (pictures 1-4) or short, wavy or tangled DNA molecules may affect the result of the test.

Do not use the product after the expiry date.

As some substances contained in this product (in low concentrations and volumes) could be harmful for health, handle the reagent with care and wear appropriate personal protective equipment. See the Material Safety Data Sheet (MSDS) for safety information.

**Reagents and materials recommended but not provided for the hybridization and the detection of the “CombHeliX® FSHD probes”**

*IMPORTANT: Since the performance characteristics of the “CombHeliX® FSHD probes” have been evaluated and validated using the reagents and materials listed below, we recommend the use of these referenced reagents and materials for an optimal diagnostic result.*

The “CombHeliX® FSHD probes” were validated with the following accompanying reagents and materials for their detection as green, red, blue and magenta fluorescent signals:

- Green fluorescent signal:
  - o 1mg/ml anti-fluorescein/Oregon Green® antibody, rabbit IgG fraction (Invitrogen, Ref: A-889). Ready-to-use.
  - o 2 mg/mL Alexa Fluor® 488 goat anti-rabbit IgG (H+L) (Invitrogen, Ref: A-11034). Ready-to-use.
- Red fluorescent signal:
  - o 1 mg streptavidin, Alexa Fluor® 594 conjugate (Invitrogen, Ref: S-11227). To be reconstituted (see reagent preparation).
  - o 0.5 mg biotinylated Anti-Streptavidin (Vector Laboratories, Ref: BA-500). To be reconstituted (see reagent preparation).
- Blue fluorescent signal:
  - o 0.5 mg aminomethylcoumarin (AMCA)-IgG Fraction Monoclonal Mouse Anti-Digoxin (Jackson ImmunoResearch, Ref: 415-155-156). To be reconstituted (see reagent preparation).

- 1 mg AMCA-AffiniPure Rat Anti-Mouse IgG (H+L) (Jackson Immunoresearch, Ref: 415-155-166). To be reconstituted (see reagent preparation).
- 2 mg/mL Alexa Fluor® 350 Goat Anti-Rat IgG (H+L) (Invitrogen, Ref: A-21093). Ready-to-use.

Note : Magenta fluorescent signal is the result of the double labeling of the corresponding probes with digoxigenin and biotin that allows the detection of a co-localized blue and red signal with the above-mentioned reagents.

- “CombHeliX® FSHD Software” (Contact Genomic Vision for further additional information).

Other reagents and materials are:

- Autoclaved distilled water
- Deionized formamide
- 20X SSC
- 70%, 90% and 100% ethanol
- BlockAid™ blocking solution
- Tween® 20
- 1X PBS
- Variable micropipette (1 µL - 200 µL)
- Tweezers
- Microscope slide
- Co-denaturation and hybridization instrument (e.g. Hybridizer, Dako)
- Ceramic coverslip tray
- Humidified chamber
- Incubator at 37°C
- 250 mL-beaker
- ImageXpress® Micro System (Molecular Devices)

**Preparation of the reagents for the hybridization and the detection of the “CombHeliX® FSHD probes”**

*IMPORTANT: Use autoclaved distilled water for preparation of all stock and working solution.*

- Hybridization Washing Buffer (2X SSC solution): mix 100 mL of 20X SSC and 900 mL of distilled water. Store at RT.

*Note: The buffer must be warmed at 60°C prior to use.*

- Detection Washing Buffer (2X SSC/1% Tween): Mix 100 mL of 20X SSC and 10 mL of Tween® 20 with 890 mL of distilled water. To be prepared extemporaneously.
- Reconstitute the lyophilized reagents with autoclaved distilled water as follows:

	Quantity (mg)	Volume of distilled water to be added (mL)	Final concentration (mg/mL)
Biotinylated Anti-Streptavidin	0,5	1	0,5
Streptavidin, Alexa Fluor® 594 conjugate	1	1	1
AMCA-IgG Fraction Monoclonal Mouse Anti-Digoxin	0.5	0,35	1,43
AMCA-AffiniPure Rat Anti-Mouse IgG	1	0,75	1,33

- Prepare extemporaneously the mix of the different layers of detecting solutions as follows (for one test):

	Volume (µL)
<b>First layer:</b>	
- 1 mg/mL Streptavidin, Alexa Fluor® 594 conjugate	0,8
- 1 mg/ml anti-fluorescein/Oregon Green® antibody, rabbit IgG fraction	0,8
- 1,43 mg/mL AMCA-IgG Fraction Monoclonal Mouse Anti-Digoxin	0,8
- BlockAid™ blocking solution	17,6
<b>Second Layer:</b>	
- 0,5 mg/mL Biotinylated Anti-Streptavidin	0,8
- 2 mg/mL Alexa Fluor® 488 goat anti-rabbit IgG (H+L)	0,8
- 1,33 mg/mL AMCA-AffiniPure Rat Anti-Mouse IgG	0,8
- BlockAid™ blocking solution	17,6
<b>Third layer :</b>	
- 1 mg/mL Streptavidin, Alexa Fluor® 594 conjugate	0,8
- 2 mg/mL Alexa Fluor® 350 Goat Anti-Rat IgG (H+L)	0,8
- BlockAid™ blocking solution	18,4

## Procedure

### **Hybridization of the CombHeliX® FSHD probes on combed DNA**

*IMPORTANT: all the steps where formamide is used should be performed under a fume hood.*

1. Thaw the vial containing the “CombHeliX® FSHD probes” and the coverslips with combed DNA on a ceramic coverslip tray for 10 min at room temperature.
2. Dehydrate the combed DNA coverslips by dipping the ceramic coverslip tray 3 min at room temperature in successive baths of 70%, 90% and 100% ethanol.
3. Air dry the coverslips at room temperature for 10 min protecting from light.
4. For one test, transfer 10 µL of the “CombHeliX® FSHD Probes” in a new microcentrifuge tube and add 10 µL of deionized formamide to the “CombHeliX® FSHD probes”
5. Mix well and incubate at 37°C for 30 min

6. Pipette 20  $\mu$ L of the probes/formamide mix on a microscope slide.
7. Avoiding trapped bubbles, set the coverslip on the drop of hybridization solution.  
*Note: Carefully indicate the side of the coverslip in contact with the hybridization solution to avoid confusion in the subsequent steps.*  
*Note: Lay down the coverslip carefully using tweezers to avoid bubbles. Do not push and slide the coverslip once mounted, it will cause scratches on combed coverslips. Adjust the position of coverslip by gentle touching the corner if needed.*  
*Note: Protect the hybridized slides from light from this step.*
8. Co-denature the combed DNA on the coverslip and the probes for 5 min at 90°C in the humidified chamber of a hands-free co-denaturation and hybridization instrument.
9. Incubate for 16-20 hrs at 37°C in the humidified chamber of the co-denaturation and hybridization instrument.
10. Remove the coverslip from the microscope slide and place on a ceramic coverslip tray in a 250 mL-beaker containing the pre-warmed (60°C) Hybridization Washing Buffer.  
*Note: if the coverslip is stuck to the slide, add a drop of Hybridization Washing Buffer 1 around the coverslip and wait until the coverslips floats. When removing coverslip, slide the coverslips gently until one of corner is out on slide to avoid scratches, then gently “peel off” the coverslip with tweezers.*
11. Wash the coverslip three times in pre-warmed (60°C) Hybridization Washing Buffer for 5 min each at 60°C.  
*Note: Do not let the coverslip dry. It must be totally immersed in the solution.*

#### **Detection of the CombHelix® FSHD probes**

1. Pipette 20  $\mu$ L of the first layer of detecting solution on a microscope slide and set the coverslip on the drop of detecting solution.  
*Note: Make sure to place the hybridized side of the coverslip in contact with the detecting solution.*
2. Incubate the slide for 20 min in a humidified chamber at 37 °C.
3. Remove the coverslip from the microscope slide and place it on a ceramic coverslip tray in a 250 mL-beaker containing the Detection Washing Buffer.
4. Wash the coverslip three times in the Detection Washing Buffer for 3 min each at room temperature with gentle agitation.
5. Pipette 20  $\mu$ L of the second layer of detecting solution on a microscope slide and set the coverslip on the drop of detecting solution.

*Note: Make sure to place the hybridized side of the coverslip in contact with the detecting solution.*

6. Repeat steps 2 to 4
7. Pipette 20  $\mu$ L of the third layer of detecting solution on a microscope slide and set the coverslip on the drop of detecting solution.

*Note: Make sure to place the hybridized side of the coverslip in contact with the detecting solution.*

8. Repeat steps 2 to 4
9. Wash the coverslip in 1xPBS for 3 min at room temperature with gentle agitation.
10. Dehydrate the coverslip by dipping the ceramic coverslip tray 3 min in successive baths of 70%, 90% and 100% ethanol.
11. Air dry the coverslip at room temperature for 10 min protecting from light.
12. The coverslip can be stored at 4°C protecting from light until observation.

### **Visualization of the fluorescent signals**


For the *in vitro* diagnostic of FSHD, the hybridized coverslips must be scanned with the ImageXpress® Micro System, an automated image scanning system from Molecular Devices, equipped with a 40X objective and the following filters:

<b>Fluorophore</b>	<b>Excitation</b>	<b>Emission</b>
Green	482 $\pm$ 35 nm	536 $\pm$ 40 nm
Red	562 $\pm$ 40 nm	624 $\pm$ 60 nm
Blue	377 $\pm$ 50 nm	447 $\pm$ 60 nm

The hybridization signals appear as multicolor fluorescent signal arrays. The 3.3kb-D4Z4 repeats appear in green while signals, which distinguish 4q and 10q chromosomes and qA and qB haplotypes, appear in red, blue and magenta with a specific size, spacing and order (as indicated in the figure below). The different fluorescent signals that might be detected are as follows:

- Chromosome 4q35.2 is identified by a centromeric 18kb-red signal, a -1.5 kb Dux4 green signal, a 10kb-magenta signal, a 5kb space and a 20kb-blue signal;
- Chromosome 10q26.1 is identified by a centromeric 15kb-blue, a 11kb space, a 10-kb magenta, a 5kb space and a 20kb-blue signal;
- Haplotype qA appears as a 2 telomeric 6kb red signals
- Haplotype qB appears as a 2 telomeric 6kb blue and red signals. An additional genomic region of unknown chromosomal localization (Chr Z), which is apparently unrelated to FSHD, is also detected as 2 blue fluorescent signals surrounding by a red and a magenta signal.

For the diagnostic of FSHD, the detection, the classification and the measurement of the different fluorescent signals (especially the length of D4Z4 repeat array) have to be performed using the “CombHeliX® FSHD software”. Please refer to the Instructions for Use of the “CombHeliX® FSHD software” or contact the manufacturer for additional information.

### Interference

The “CombHeliX® FSHD probes” generates an additional signal on Chr 3 (indicated as Chr 3 in the figure above) but does not interfere with the results since its specific fluorescent pattern allows its identification. These signals are not included in the analysis and in the diagnostic result.

### Recommendations and limitations

The “CombHeliX® FSHD probes” are designed for *in vitro diagnostic use* (according to the EU directive 98/79/EC) and have to be analyzed with the “CombHeliX® FSHD Software”. Interpretation of results must be made within the context of the patient’s clinical history by a qualified pathologist. A medical decision cannot be made based on the results of this assay alone. The “CombHeliX® FSHD probes” cannot be used for pre-natal diagnostic or for other diagnostic use (e.g. other muscular disorders). Genomic Vision does not assume any responsibility for improper application of this product.

### Training

Training is mandatory for performing this assay. Genomic Vision will provide training in specimen preparation, assay procedure, and interpretation of Molecular Combing testing

of FSHD for inexperienced users. It is also recommended that a laboratory that has previously received training but now has new personnel performing the assay request training for the new users.

### **Performance characteristics**

The performance characteristics have been established on 69 patients collected and analyzed by two independent investigation centers at l'Hôpital La Timone (Marseille, France) and at the Institute of Human Genetics (Wurzburg, Germany).

The manufacturer ensures that this product meets the performance characteristics established on the intended sample:

- High accuracy of the number of D4Z4 repeats:
  - $\pm 1$  D4Z4 repeat = 89 % ;
  - $\pm 2$  D4Z4 repeats = 98 %.
- High sensitivity: >99%. The detection of a contracted 4qA allele (<10 D4Z4 copies) leads to a FSHD1 positive diagnostic.
- High specificity: >99%. Absence of a contracted 4qA allele (>10 D4Z4 copies) results to a FSHD1 negative diagnostic.

### **References**

1. van Deutekom, J.C., et al., Hum Mol Genet, 1993. **2**(12): p. 2037-42.
2. Wijmenga, C., et al., Nat Genet, 1992. **2**(1): p. 26-30.
3. Winokur, S.T., et al., Chromosome Res, 1994. **2**(3): p. 225-34.
4. Hewitt, J.E., et al., Hum Mol Genet, 1994. **3**(8): p. 1287-95.
5. Lemmers, R.J., et al., Nat Genet, 2002. **32**(2): p. 235-6.
6. Lemmers, R.J., et al., Am J Hum Genet, 2004. **75**(6): p. 1124-30.
7. van Geel, M., et al., Genomics, 2002. **79**(2): p. 210-7.
8. Nguyen, K., et al., Ann Neurol, 2011. **70**(4): p. 627-33.
9. Bensimon, A., et al., Science, 1994. **265**(5181): p. 2096-8.
10. Lebofsky, R. and A. Bensimon. Brief Funct Genomic Proteomic, 2003. **1**(4): p. 385-96.

## Ordering information and related products

Product	Content	Reference
CombHelix® FSHD Hybridization kit	10 hybridizations	FSH-HYB-001
CombHelix® FSHD Software	1 permanent license	FSH-ASW-001
CombHelix® DNA Extraction kit *	100 extractions	EXT-001
Molecular Combing System *	Combing device	MCS-001
Teflon Reservoir (for DNA solutions)	1 piece, reusable	TEF-001
Coverslip holder	1 piece, for 2 coverslips	CLI-001
Disposable Reservoirs *	Pack of 10 pieces	RES-001
Reservoir Supports for MCS	Pack of 2 pieces	SUP-001
Reservoir Bench Holder	1 piece, 10 positions	POR-001
Silanised Coverslips *	Box of 50 pieces	COV-001

\* Products for general laboratory use that have been validated for Molecular Combing but none specific validation for *in vitro* diagnostic use has been conducted.

Our experts are available to answer your questions

<http://www.genomicvision.com>

[support@genomicvision.com](mailto:support@genomicvision.com)

### Patent

This product or the use of this product is subject to proprietary rights (EP2007/059299, IB2009/007197). The Molecular Combing technology and products are covered by patents (FR2716206, FR 2716263, FR 2737574, FR 2755149) owned by Genomic Vision S.A.

### **Trademarks:**

CombHelix is a registered trademark of Genomic Vision S.A.

BlockAid is a trademark of Molecular Probes, Inc.

Alexa Fluor and Oregon Green are registered trademarks of Molecular Probes, Inc.

ImageXpress Micro System is a registered trademark of Molecular Devices, LLC.

Tween is a registered trademark of Croda International PLC

Genomic Vision S.A.  
80-84 rue des meuniers  
92220 Bagneux  
France  
[www.genomicvision.com](http://www.genomicvision.com)

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CombHeliX est une marque déposée par Genomic Vision S.A.
BlockAid est une marque déposée par Molecular Probes, Inc.
Alexa Fluor et Oregon Green sont des marques déposées par Molecular Probes, Inc.
ImageXpress Micro System est une marque déposée par Molecular Devices, LLC.